oobrush®

INSTRUCTIONS - BATTERY-POWERED TOOTHBRUSH

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**DENTAL
REVOLUTION**

USING THE TOOTHBRUSH

1 Insert the batteries as shown in **Figure A**. Operates with two 1.5V AA (LR6) batteries. For improved performance, we recommend that you use alkaline batteries.

2 For your comfort and to avoid spraying, the toothbrush features a continuous pressure on-off switch. To turn the toothbrush on, maintain continuous pressure on the on-off switch.

CAUTION: Keep the brush away from hair to avoid any risk of entanglement.

3 Dampen the brush head, apply a small amount of toothpaste between the two brushes and switch the unit on for a few seconds to spread the toothpaste evenly. Place the unit in your mouth so that the two brushes are on opposite sides of your teeth and move slowly along the dental arch with a slight back and forth movement (**Figure B**). Thus, you brush simultaneously the inside and outside surfaces of the teeth while always respecting the correct brushing method from gum to tooth.

OObrush massages your gums thoroughly, which reinforces their vitality. When using the toothbrush for the first time, you may experience a temporary feeling of gum irritation or tickling.

This is a normal phenomenon associated with the massaging action of the gums and should disappear naturally.

4 There are several ways to use your toothbrush:

Figure C: Usual brushing position: simultaneous brushing of the inside and outside surfaces of the teeth.

Figure D: Brushing of top surfaces of the teeth, in particular for the molars.

Figure E: Brushing the outside surfaces of both the top and bottom teeth with teeth clenched.

5 To clean the brush, turn the appliance on and rinse the brush head thoroughly under running water for a few seconds.

6 We recommend you replace the brush head of your toothbrush every three months, as shown (**Figure F**).

7 Dispose of used batteries according to the applicable regulations. At the end of the product's life, please bring it to an approved recycling center or to your retailer.

**GUARANTEE**

The guarantee is valid two years starting on the date of purchase. The guarantee

stipulates that the product will be replaced free of charge in the case of design or manufacturing defects.

The guarantee does not cover products that have been opened by unauthorized persons. It does not cover normal battery use and ordinary wear. It does not apply to damages resulting from transport, dropping, or improper handling of the product. It does not cover inappropriate uses of the product.

For all claims during the guarantee period, return the unit (charger and brush) along with the receipt to the dealer.

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