# oobrush®

INSTRUCTIONS - RECHARGEABLE TOOTHBRUSH

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[www.oobrush.com](http://www.oobrush.com)

**1. IMPORTANT SAFEGUARDS**

When using electrical products, especially when children are present, certain basic safety precautions should be followed. Please read the following instructions before using your power toothbrush..

**WARNING**

To reduce the risks of burns, electrocution, fire or injury, please follow the instructions below:

1. The supply cord cannot be replaced. If the cord is damaged the appliance should be scrapped.
2. Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision.
3. Children shall not play with the appliance.
4. Do not place or store charger where it can fall or be pulled into a tub or sink.
5. Do not submerge charger in water or any other liquid.
6. Do not touch the charger if it has been in contact with water. Unplug it immediately.
7. Do not use while bathing..
8. Exercise close supervision when the unit is used by children or disabled persons.
9. Use this product only for its intended use as described in this manual..
10. Do not use attachments not recommended by the manufacturer.
11. Never use this product if it has a damaged cord or plug, if it is not working properly, if it has fallen or been damaged, or if the charger has come in contact with water. The power cord cannot be replaced; if damaged, it must not be used.
12. Keep the cord away from heated surfaces.
13. Never drop or insert foreign objects into any opening.

**SAVE THESE INSTRUCTIONS FOR THE PRODUCT'S ENTIRE LIFE**

## 2. GETTING STARTED AND CHARGING

The charger adapts automatically to the different voltages and frequencies of all countries (110V to 240V and 50Hz to 60Hz).

The toothbrush must be fully charged before first use. To charge the brush, place it in the charger and plug the charger into an electrical outlet. The charge indicator light goes on. The initial charge takes at least 24 hours.

For daily use, the toothbrush can be stored on the charger without any risk of overcharging. Once fully charged, the unit will run 20 minutes on a single charge.

To maintain the maximum capacity of the battery, unplug the charger and discharge the brush by regular use at least every six months.

## 3. USING THE TOOTHBRUSH

For your comfort and to avoid spraying, the toothbrush features a continuous pressure on-off switch. To turn the toothbrush on, maintain continuous pressure on the on-off switch.

**CAUTION:** Keep the brush away from hair to avoid any risk of entanglement.

Dampen the brush head, apply a small amount of toothpaste between the two brushes and switch the unit on for a few seconds to spread the toothpaste evenly. Place the unit in your mouth so that the two brushes are on opposite sides of your teeth and move slowly along the dental arch with a slight back and forth movement (**Figure A**). Thus, you brush simultaneously the inside and outside surfaces of the teeth while always respecting the correct brushing method from gum to tooth.

OObrush massages your gums thoroughly, which reinforces their vitality. When using the toothbrush for the first time, you may experience a temporary feeling of gum irritation or tickling. This is a normal phenomenon associated with the massaging action of the gums and should disappear naturally.

**There are several ways to use your toothbrush:**

**Figure B:** The usual brushing position brushes simultaneously the inside and outside surfaces of the teeth.

**Figure C:** The 45° angle position brushes the top surfaces of the teeth, in particular for the molars.

**Figure D:** The 90° angle position brushes the outside surfaces of both the top and bottom teeth, while keeping teeth clenched.

To clean the brush, put the unit on and rinse carefully the brush head under running water during several seconds (**Figure E**).

## 4. REPLACING THE BRUSH HEAD

We recommend you replace the brush head of your toothbrush every three months, as shown in **Figure F**.

## 5. RESPECTING THE ENVIRONMENT

The brush uses a non-replaceable rechargeable battery. Follow the regulations in force for the disposal of the brush at the end of its useful life. Dispose of it at an appropriate collection site so that it can be disposed of in an environmentally-safe manner. To remove the battery, open the body of the brush and cut the power leads. To avoid any danger, make sure that the battery is fully discharged before removal.

## 6. GUARANTEE

The guarantee is valid two years starting on the date of purchase. The guarantee stipulates that the product will be replaced free of charge in the case of design or manufacturing defects.

The guarantee does not cover products that have been opened by unauthorized persons. It does not cover ordinary wear and does not apply to damages resulting from transport, dropping, or improper handling of the product. It does not cover inappropriate uses of the product.

For all claims during the guarantee period, return the unit (charger and brush) along with the receipt to the dealer.

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